



EBERLE

Shaver System C3



EBERLE

Made in Germany

Shaver System C3



Control unit


- Panel to control motor-powered EBERLE C3 handpieces
- Resistive 5-inch touch display with automatic brightness adjustment
- Additional softkeys to set basic functions
- Selectable start-up behaviour for handpieces (soft, medium and hard)
- Freely programmable memory areas to define user-related or use-oriented settings
- User-specific assignment of control elements on handpiece and footswitch
- Rotational speed range 300 – 10.000 min-1
- Tapping-mode for Universal handpiece

Accessories



- Handpieces of various housing materials, with and without manual operation
- Universal handpiece with various attachments for drilling and sawing
- Single-use and reusable shaver blades
- 5-button footswitch

Video:







Control unit C3





Universal handpiece





○ **Shaver System**

..... **Control unit**
 article no. 913055

..... **Footswitch**
 article no. 913057

..... **Handpiece**
 article no.
 913066 8k
 913070 10k
 320 g + 4 m cable
 casing material aluminum
 article no.
 913068 8k
 913080 10k
 450 g + 4 m cable
 casing material stainless steel

Handpiece, manual operation
 article no.
 913067 8k
 913071 10k
 380 g + 4 m cable
 casing material aluminum
 article no.
 913069 8k
 913081 10k
 540 g + 4 m cable
 casing material stainless steel

Small Joint Handpiece
 article no.
 913065SJ
 185 g + 4 m cable
 casing material aluminum

Control unit

- Panel to control motor-powered EBERLE C3 handpieces
- 5-inch touch display
- Automatic brightness adjustment
- Softkeys for basic functions
- Adjustable oscillation
- Oscillation mode: revolutions or frequency
- User-specific assignment of the control elements for the footswitch as well as for the handpiece
- Individual acceleration behaviour
- User name editable
- Ergonomic 5-button footswitch

Handpiece

- ergonomically shaped handpiece, autoclavable
- micro-processor controlled drive motor

Handpiece, manual operation

- ergonomically shaped handpiece, autoclavable
- micro-processor controlled drive motor
- manual operation of several functions



Speed ranges	Handpiece 8k	Handpiece 10k
Cutting	300–4000 min ⁻¹	500–4000 min ⁻¹
Oscillation	500–4000 min ⁻¹	500–4000 min ⁻¹
Burring	500–8000 min ⁻¹	500–10000 min ⁻¹

Small Joint Handpiece

article no. 913065SJ

- ergonomically shaped handpiece, autoclavable
- micro-processor controlled drive motor
- speed ranges:
 - rotation 500 – 6000 min⁻¹
 - oscillation 500 – 4000 min⁻¹



Footswitch

article no. 913057

- user-friendly, ergonomic 5-button footswitch



Cleaning adapter

for EBERLE shaver handpieces


article no. 100020



Brushes

for EBERLE handpieces

Cleaning brush, bent-over tip
length: 150 mm
article no. 011001-301



Cleaning brush with eyelet
length: 150 mm
article no. 011001-302



Cleaning brush, rounded tip
length: 300 mm
article no. 011001-300





**Universal handpiece
for manual operation**

- casing material aluminium
- weight 750 g + 4 m cable
- ergonomically shaped handpiece, autoclavable
- micro-processor controlled drive motor
- manual operation of anticlockwise/clockwise rotation
- speed range: 0 – 1500 min⁻¹, infinitely variable

○ **Universal handpiece**
article no. 913075

..... **Wire driver adapter**
article no. 913026
article no. 913034

..... **AO Drill adapter**
article no. 913044

..... **Drill chuck adapter**
article no. 913027
article no. 913019
article no. 913022

..... **Sagittal saw adapter**
article no. 913024

**Cleaning adapter
for Universal handpiece**

article no. 100025



**Modular attachments
for Universal handpiece**

Wire driver adapter

article no. 913026

- receptacle for 1.8 – 3.2 mm

article no. 913034

- receptacle for 0.6 – 1.8 mm



AO Drill adapter

article no. 913044

- receptacle for standard AO



Drill chuck adapter

article no. 913027

- receptacle for 0.7 – 7.4 mm



Drill chuck adapter, high-power

article no. 913019

- receptacle for 0.7 – 7.4 mm
- rotational speed 0 – 500 min⁻¹



Drill chuck adapter

article no. 913022

- receptacle for 0 - 3 mm, keyless



Sagittal saw adapter

article no. 913024

- oscillation frequency 0 – 22.500 cycl./min.



**Cleaning adapter
for modular attachments**

article no. 100026













Variety of the Shaver Blades










Shaver Blades reusable






effective length 130 mm
Packing unit 5 pieces

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Full Radius Resector	9131135	9131142	9131148	9131155	
 Full Radius Resector Sharp Cut		9132742			
 Aggressive Full Radius Resector	9131335	9131342	9131348	9131355	
 Aggressive Full Radius Resector Plus		9133442			
 Double Full Radius Resector		9132542		9132555	
 Double Aggressive Full Radius Resector		9132342		9132355	
 Conical Double Full Radius Resector		9132442			
 Full Radius Resector Small Window		9132242		9132255	
 Meniscus Cutter	9131235	9131242		9131255	
 Aggressive Meniscus Cutter	9131535	9131542	9131548	9131555	

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Aggressive Meniscus Cutter Plus		9132142			
 Double Meniscus Cutter		9132642		9132655	
 Aggressive Meniscus Side Cutter	9131435	9131442		9131455	
 Round Burr		9131042 Ø burr 3.0 mm		9131055 Ø burr 4.0 mm	9131065 Ø burr 5.0 mm
 Oval Burr				9131955 Ø burr 4.0 mm	9131965 Ø burr 5.0 mm
 Aggressive Oval Burr				9133855 Ø burr 4.0 mm	9133865 Ø burr 5.0 mm
 Tapered Burr				9131855 Ø burr 4.0 mm	9131865 Ø burr 5.0 mm











effective length 130 mm, 15° bent

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Full Radius Resector cutting window inward		9133042		9133055	
 Aggressive Full Radius Resector cutting window inward		9133242		9133255	
 Aggressive Meniscus Side Cutter cutting window outward		9133142		9133155	

Shaver Blades reusable




for hip arthroscopy, effective length 210 mm
Packing unit 5 pieces

		Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
	Full Radius Resector Sharp Cut tubular		9132742H			
	Aggressive Full Radius Resector tubular		9131342H	9131348H	9131355H	
	Meniscus Cutter tubular		9131242H		9131255H	
	Aggressive Meniscus Cutter tubular		9131542H	9131548H	9131555H	
	Aggressive Meniscus Cutter Plus tubular		9132142H			
	Round Burr				9131055H Ø burr 4.0 mm	
	Oval Burr				9131955H Ø burr 4.0 mm	
	Aggressive Oval Burr				9133855H Ø burr 4.0 mm	







for hip arthroscopy, effective length 210 mm, 15° bent

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Aggressive Full Radius Resector cutting window inwards		9133242H		9133255H	




for small joint surgery, effective length 70 mm

	Ø 2.0 mm article no.	Ø 2.5 mm article no.	Ø 3.5 mm article no.
 Full Radius Resector	9131120SJ	9131125SJ	9131135SJ
 Aggressive Full Radius Resector	9131320SJ	9131325SJ	9131335SJ
 Round Burr		9131025SJ Ø burr 1.8 mm	9131035SJ Ø burr 2.5 mm



Cleaning adapter

for reusable
C3 compatible shaver blades
article no. 100001













Cleaning set








article no. 900019

Shaver Blades single-use






effective length 130 mm, sterile packaging
Packing unit 5 pieces

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Full Radius Resector	9141135	9141142	9141148	9141155	
 Full Radius Resector Sharp Cut		9142742		9142755	
				Burr speed	
 Aggressive Full Radius Resector	9141335	9141342	9141348	9141355	
 Aggressive Full Radius Resector Plus		9143442			
 Double Full Radius Resector		9142542		9142555	
 Double Aggressive Full Radius Resector		9142342		9142355	
 Conical Double Full Radius Resector		9142442			
 Full Radius Resector Small Window		9142242		9142255	
 Meniscus Cutter	9141235	9141242		9141255	
 Aggressive Meniscus Cutter	9141535	9141542	9141548	9141555	

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Aggressive Meniscus Cutter Plus		9142142			
 Double Meniscus Cutter		9142642		9142655	
 Round Burr		9141042 Ø burr 3.0 mm		9141055 Ø burr 4.0 mm	9141065 Ø burr 5.0 mm
 Oval Burr				9141955 Ø burr 4.0 mm	9141965 Ø burr 5.0 mm
 Aggressive Oval Burr				9143855 Ø burr 4.0 mm	9143865 Ø burr 5.0 mm
 Tapered Burr				9141855 Ø burr 4.0 mm	9141865 Ø burr 5.0 mm
 Crown Cutter				9143755	





effective length 130 mm, 15° bent, sterile packaging

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Full Radius Resector cutting window inward		9143042		9143055	
 Aggressive Full Radius Resector cutting window inward		9143242		9143255	
 Aggressive Meniscus Side Cutter cutting window outward		9143142		9143155	




for hip arthroscopy, effective length 210 mm, sterile packaging
Packing unit 5 pieces

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Full Radius Resector Sharp Cut tubular		9142742H			
 Aggressive Full Radius Resector tubular		9141342H	9141348H	9141355H	
 Meniscus Cutter tubular		9141242H		9141255H	
 Aggressive Meniscus Cutter tubular		9141542H	9141548H	9141555H	
 Aggressive Meniscus Cutter Plus tubular		9142142H			
 Round Burr				9141055H Ø burr 4.0 mm	
 Oval Burr				9141955H Ø burr 4.0 mm	
 Aggressive Oval Burr				9143855H Ø burr 4.0 mm	




for hip arthroscopy, effective length 210 mm, 15° bent, sterile packaging

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Aggressive Full Radius Resector cutting window inwards		9143242H		9143255H	






for hip arthroscopy, effective length 210 mm, 20° bent, sterile packaging

	Ø 3.5 mm article no.	Ø 4.2 mm article no.	Ø 4.8 mm article no.	Ø 5.5 mm article no.	Ø 6.5 mm article no.
 Round Burr		9143642H			
		Ø burr 5 mm			



for small joint surgery, effective length 70 mm, sterile packaging

	Ø 2.0 mm article no.	Ø 2.5 mm article no.	Ø 3.5 mm article no.
 Full Radius Resector	9141120SJ	9141125SJ	9141135SJ
 Aggressive Full Radius Resector	9141320SJ	9141325SJ	9141335SJ
 Round Burr		9141025SJ	9141035SJ
		Ø burr 1.8 mm	Ø burr 2.5 mm

EBERLE GmbH & Co. KG

Dachsteinstraße 2
75449 Wurmberg
Germany

Phone +49 (0) 7044 9611-288

sales@eberle-med.com
www.eberle-med.com

